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Leybold Handblätter Physik, Teil 1 (597 621)

The spectrometer and goniometer (467 23) is used

- to measure emission and absorption spectra, and
- to determine characteristic data such as refracting angle, angle of minimum deviation, refractive index and mean dispersion.

The device is supplied with a prism; however, when equipped with a grating (e.g. line grating, 471 51) this apparatus can also be used as a grating spectrometer.

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1 Safety notes

- This is a precision optical instrument: protect it from mechanical shocks!
- Max. voltage for lamp of illumination device: 8 V!


Fig. 1


Fig. 2

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2 Scope of supply, description, technical data

2.1 Scope of supply

- Ⓐ Spectrometer base assembly
- Ⓑ Telescope
- Ⓒ Slit tube
- Ⓓ Illumination device
- Ⓔ Table for Ⓕ and Ⓖ
- Ⓕ Holder with flat glass plate (or for mounting a grating, e.g. 471 51)
- Ⓖ Flint glass prism on holder

2.2 Description, technical data (see Fig. 2)

- ① Telescope
 - (1.1) Adjusting screw for compensating alignment errors between telescope ① and slit tube ㉑
- ② Focus screw
- ③ Sliding Gauss eyepiece
 - (3.1) Aperture for light from ④
- ④ Illumination device
 - (4.1) Arresting screw for illumination device
- ⑤ Adjusting screws for lateral displacement of telescope ①
- ⑥ Height-adjustment screw for telescope ①, lockable
- ⑦ Levelling screw for prism table ㉒
- ⑧ Arresting screw for prism table ㉒
- ⑨ Arresting screw for disk scale ⑩
- ⑩ Disk scale
- ⑪ Fine adjustment for telescope rotation
- ⑫ Arresting screw for telescope ① (covered in drawing)
- ⑬ Verniers
- ⑭ Vernier magnifiers
- ⑮ Height-adjustment screw for slit tube
- ⑯ Adjusting screw for lateral displacement of slit tube ㉑
- ⑰ Arresting screw for slit width
- ⑱ Adjustable slit stop
- ⑲ Adjustable slit
- ⑳ Micrometer screw for slit width
- ㉑ Slit tube
 - (21.1) Adjusting screw for compensating alignment errors with respect to telescope
- ㉒ Table for ㉓ or ㉔
- ㉓ Holder with flat glass plate (adjusting aid) or for mounting a grating (grating spectrometer)
- ㉔ Holder with prism (prism spectrometer)

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Technical data

Telescope:	focal length $f = 160$ mm; aperture 16 mm; Gauss eyepiece with illuminated cross wires; focused with knurled screw; freely pivoting with respect to axis of apparatus
Slit tube:	focal length $f = 160$ mm; aperture 16 mm; adjustable slit width
Slit:	width adjustable via micrometer screw.
Disk scale:	diameter 180 mm full and half-degree scale divisions; freely turning with reference to axis of apparatus
Vernier:	two verniers for reading disk scale with divisions 0 - 30, opposite each other; with reading magnifiers
Prism:	flint glass prism, on holder; base: 33 mm; height 22 mm Refractive index n_D : 1.620 mean dispersion $n_F - n_D$: 0.017 resolution $\lambda/\Delta\lambda$: approx. 3300
Illumination device for cross wires:	with lamp 8 V/0,15 A (additionally required: voltage source e.g. transformer 6 V/12 V, 562 73).
Adjusting aid:	flat glass plate on holder

3 Operation

3.1 Initial assembly

Mount the telescope, illumination device, slit tube and prism table on the spectrometer base assembly as shown in Fig. 2.

3.2 Adjusting the apparatus

In order to conduct accurate measurements, you need to adjust the apparatus carefully.

The slit and cross wires of the telescope must be in the focal planes of the corresponding objective lenses (telescopic beam path).

Align the slit and the prism surfaces parallel to the rotary axis of the telescope.

Some of the working steps for adjusting the apparatus and measuring the line spectra are easier to carry out when the room is darkened somewhat.

3.2.1 Preliminary adjustment

The inclination of the prism table can only be varied to a limited extent. In order to ensure sufficient play for fine adjustment, align the prism table by eye as nearly horizontal as possible before beginning with the preliminary adjustment.

- Align the telescope ①, prism table ② and slit tube ③ horizontally by eye.
- Align the telescope and slit tube with the middle using the adjusting screws ⑤, ⑥ for lateral displacement and lock them in position. Do not loosen the adjusting screws too much on any one side, as these support the telescope and the slit tube.

3.2.2 Adjusting the telescope to infinity

Note: although experimenters with imperfect vision can see distant objects clearly with the telescope, by definition the telescope is not exactly focused on infinity in this case. However, it is still possible to carry out precise measurements when the telescope is (approximately) focused on infinity and the same experimenter adjusts the slit tube. Other experimenters who wish to view the spectra may only focus the setup by sliding the eyepiece ③.

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- Remove eyepiece ③, attach illumination device ④ to telescope and remount the eyepiece with the light aperture (3.1) facing downward.
- Focus the cross wires by sliding the eyepiece ③ in the eyepiece tube and align them if necessary. Make sure that the opening for the illumination device on the eyepiece is still facing downward.
- Focus the horizontally aligned telescope on a distant object (500 m) with the focusing screw ②.

When the apparatus is correctly adjusted, the image of the observed object and the cross wires will be in the focal plane of the objective; thus there should be no parallax between the object and the cross wires.

3.2.3 Adjusting der illumination device

- Align the telescope with the slit tube (slit opened slightly).
- Connect the illumination device ④ to 6 V to max. 8 V
- Use the arresting screw (4.1) to mount the illumination device on the telescope so that the inside of the slit is well illuminated, without disturbing the eyepiece.


Fig. 4

3.2.4 Adjusting the optical axis of the telescope perpendicular to the spectrometer axis

- Place the flat glass plate on holder ㉑ in the center of the prism table at an angle of 45° to slit tube ㉑ so that the imaginary line between two leveling screws of the prism table is parallel to the surfaces of the flat glass plate (see Fig. 3).
- Align telescope ① perpendicular to one face of the flat glass plate so that the cross wires are reflected in the surface.
- Place the horizontal cross wire exactly over its mirror image; correct half the distance using the height adjustment screw ⑥ of the telescope and the other half with the leveling screw ⑦ of the prism table indicated in Fig. 3.
- Repeat the following two steps until the horizontal cross wire and its mirror image remain coincident on both sides of the flat glass plate:
 - 1) Pivot the telescope by 180° as shown in Fig. 3 so that the cross wires are reflected on the other side of the flat glass plate.
 - 2) Check whether the cross wire and its image match. If not, correct half the distance using the height adjustment screw of the telescope ⑥ and the other half with the leveling screw ⑦ of the prism table as described above.
- Lock height adjustment screw ⑥ of the telescope using the lock nut.

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- Remove the flat glass plate with holder from the prism table.
- Disconnect the illumination device from the voltage supply.

3.2.5 Adjusting the slit tube

- Illuminate the slit from the outside e.g. with light from a light bulb or one of the spectral lamps.
- Aim the telescope at the slit tube and open the slit somewhat with the micrometer screw for slit width ⑳.
- Set the slit with adjustable slit stop ⑱ to a suitable observable slit height.
- Use the height adjustment screw of the slit tube ⑮ to position the slit so that it is centered on the horizontal cross wire and arrest it in position.
- Loosen the arresting screw for slit width ⑰ and slide the slit tube in the direction of the arrow (see Fig. 2) so as to obtain a clear image.
- By turning the tube, align the slit vertically so that it is parallel to the vertical cross wire, and retighten the arresting screw for slit width ⑰.

3.2.6 Aligning the prism faces parallel to the axis of rotation

- Pivot the telescope to an acute angle to the slit tube (see Fig. 4) and lock it with arresting screw ⑫.
- Place the prism on holder ⑭ on the prism table as shown in Fig. 4 so that one of the prism surfaces is parallel to the imaginary line between two leveling screws of the prism table.
- Turn the prism table until the slit image reflected by a prism face is visible in the telescope, and lock with arresting screw for prism table ⑧.
- Use the back leveling screw ⑦ of the prism table (marked in Fig. 4) to position the reflected slit in the center of the cross wires.
- Repeat the following two steps until the reflection of the slit does not change vertically when the prism table is turned:
 - 1) Loosen arresting screw ⑧ for the prism table, turn the prism table until the slit image is reflected on the other prism face and lock it with arresting screw for prism table.
 - 2) Use the back leveling screw of the prism table to position the reflected slit in the center of the cross wires.

Adjustment is now complete.

3.3 Determining the refracting angle and the minimum deviation

Note:

To minimize reading errors and compensate for any possible eccentricity of the circular scale division with respect to the axis of rotation, the spectrometer is equipped with two verniers arranged opposite each other. The values of the two verniers are added together, the mean value found and added to the value shown on the disk scale.

Refracting angle

Align the prism table ㉒ as described in section 3.2.6 and turn it so that the refracting edge of the prism is facing the center of the slit tube ㉑ (Fig. 5). Arrest prism table ㉒ with adjusting screw ⑧.

Pivot telescope ① until slit ⑲ is reflected in one face of the prism and visible in the center of the cross wires. Use the two verniers ⑬ to read off the corresponding value on disk scale ⑩.

Then pivot the telescope until slit ⑲ is reflected in the other face of the prism.


Fig. 5

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
Fig. 6

Read off the corresponding value on disk scale ⑩. The difference between the two readings is double the refracting angle.

Angle of minimum deviation

Illuminate slit ⑲ with sodium light (e.g. sodium spectral lamp 451 111).

To find the angle of minimum deviation, position the prism so that the beam path is symmetrical for a medium wavelength of the spectrum (about 560 - 600 nm), and the deviation is thus at a minimum. This results in maximum spectral resolution.

Slowly turn prism table ⑳ and observe the shift of the slit image with telescope ① until you find the minimum position, i.e. until the slit image only moves to one side regardless of the direction of rotation of the prism table ⑳.

Lock the telescope ① and prism table ⑳ in the minimum position using the respective arresting screws ⑫ and ⑧ and read off the value from the disk scale ⑲.

Remove the prism from the prism table ⑳ and aim telescope ① at the slit ⑲. Lock the telescope ① with arresting screw ⑫ and read off the value from the disk scale ⑲. The difference between the two readings is the angle of minimum deviation.